


SpineJack[®] system

controlled anatomical
restoration



Product overview

SpineJack[®] system

controlled anatomical restoration



It is common practice for a fracture of any weight-bearing joint to first be reduced, then stabilised. The concept of SpineJack[®] system is to achieve a biomechanical restoration to allow early mobilisation and weight-bearing.

Anatomical reduction means restoration of the geometry of the whole vertebral body, which includes the cortical rings and endplates.

Anatomical restoration consists of achieving sagittal and coronal balance. These are key factors for kyphosis management and reduction of adjacent level fractures.^{1-5,7}

Vertebral endplate restoration has been described as having a positive influence on disc creeping, disc degeneration, compensatory curvatures or facet joint arthritis.^{1,6,8-11}

Studies have shown a correlation between vertebral deformation and clinical problems such as post-traumatic kyphosis, which has been depicted as a serious post-traumatic deformity.^{2,12}

Within this context, the SpineJack[®] system implant provides clinicians with a fully controlled and comprehensive treatment option for vertebral compression fractures (VCFs). This enables anatomical restoration and stabilisation.

Anatomical restoration

To achieve **anatomical restoration** after a VCF, consider the following:

- Controlled uni-directional craniocaudal expansion to restore **sagittal angulation**
- Adaptation of implant's expansion to restore **coronal angulation**
- Adaptation of implant's positioning for **endplate restoration**





Height restoration measurement

The anatomical restoration achieved by SpineJack® system has been shown using 3D reconstructions of pre and post-op CT scans.¹³ The superimposition of these images allows accurate measurement of vertebral body height restoration.

Colour-coded **3D mapping** provides visualisation of the amount of vertebral body height restoration.

Green: lowest restoration area
Red: highest restoration area



A continuum in control

Implant positioning

Controlled by specific instrumentation

Implant positioning in both sagittal and transverse planes can help achieve the best fit for the fracture's shape and the patient's anatomy.



Implant expansion

Controlled by millimetric implants expansion

Millimetric expansion of the implant can be maintained until the biomaterial is injected.






Optimal PMMA cement positioning and interdigitation

Controlled by PMMA cement fixed pathway and preservation of surrounding trabecular bone

Fixed pathway for the insertion of PMMA cement through the implant helps minimize the risk of posterior leakage. Preserving the surrounding trabeculae by a craniocaudal expansion allows for better interdigitation, thereby improving fixation and the bone healing process¹⁴⁻¹⁷.

Depending on the quality of the preserved trabecular structure, Stryker offers a range of **injectable PMMA cements**:

- Cohesion® Bone Cement (CM0300)
- VertaPlex™ HV Radiopaque Bone Cement Twin Pack (0406-622-000)
- VertaPlex™ Radiopaque Bone Cement Twin Pack (0406-422-000)



Surgical technique guide

The SpineJack® system is indicated for use in the anatomical reduction of vertebral compression fractures that may result from osteoporosis, trauma (VCF type A Magerl classification) and malignant lesions (myeloma or osteolytic metastasis).



SpineJack® system

Kit composition and implant dimensions

Preoperative planning strategy

Patient positioning and anaesthesia

Surgical steps

Step 1 - Vertebral body access

Step 2 - Implant site preparation

Step 3 - Implant insertion and expansion

Step 4 - Expander removal

Step 5 - PMMA cement preparation and injection

Kit composition

1 - Preparation kit

Guidewire (2)

Working cannula

Reamer with working cannula

Template

Cannula plug



SpineJack® preparation kits

KP004 Preparation kit Ø 4.2mm (unit)

KP001 Preparation kit Ø 5.0mm (unit)

KP058 Preparation kit Ø 5.8mm (unit)

2 - Expansion kits



SpineJack® expansion kits

KE004 Expansion kit Ø 4.2mm (unit)


KE001 Expansion kit Ø 5.0mm (unit)

KE058 Expansion kit Ø 5.8mm (unit)



Implant dimensions

Full control over positioning



The instrumentation has been developed to give complete sequential control throughout each stage of the procedure. During each step, you can control the positioning of the implant and adjust it to the desired position.




4.2mm



5.0mm



5.8mm



Preoperative planning strategy

Fracture mobility assessment

SpineJack® system is designed to treat mobile VCF's that may result from osteoporosis, trauma and malignant lesions. Assessment of the fracture's mobility prior to operating is recommended in order to maximize fracture reduction.

Implant positioning

The extent of fracture reduction depends largely on the positioning of the implant within the vertebral body. Therefore, it is recommended to map the optimal placement of the implant prior to surgery. The placement of 2 implants is often advised to achieve optimal anatomical restoration.

Vertebral dimensions


In order to ensure an optimal fit of SpineJack® implants, a CT scan of the vertebral body prior to surgery is required to confirm the adequacy of the vertebral dimensions. Access to the vertebral body requires a pedicle with a minimum diameter of 5mm: diameter KE004 + 0.8mm: 5mm.

Optimal implant fit chart

Pedicle diameter	Recommended kit
5.0-5.8mm	4.2mm
5.8-6.6mm	5.0mm
6.6mm and over	5.8mm

Tips and tricks

Consider vertebral body length to select the appropriate implant size.




Positioning and anaesthesia

Patient positioning

The patient is placed in a prone position to minimise loading on the fractured vertebra.

A hyper-lordotic position is recommended for lumbar fractures.



Anaesthesia

General, local or regional anaesthesia can be used depending on clinician preference and the patient's condition.

Tips and tricks

For all levels to be treated with SpineJack[®] system, pre-operatively assess the following on CT scans:

1. The inner diameter of the pedicle in order to define the largest size implant (including working cannula) which can potentially be inserted through the pedicle.
2. The inner diameter of the vertebral body in order to define the largest size implant which can be opened in the vertebral body.
3. The ideal positioning of the SpineJack[®] implant(s) for VCFs is described in the document "Positioning of the SpineJack[®]"¹⁸.

Controlled anatomical restoration

According to the preoperative planning strategy, a trocar is used under fluoroscopic control to determine the path to the vertebral body and to optimally position the implants (craniocaudal and mediolateral angles).

The entry point for the trocar tip should be inside the pedicle ring, close to its lateral wall, on the AP view. After advancing through the pedicle tunnel and reaching the posterior wall of the vertebral body, on the sagittal view, the tip of the trocar should be inside the pedicular ring, close to its medial wall, on the AP view.

Access to the vertebral body requires a pedicle with a minimum diameter of at least 5mm.

- Perform access procedure for both implants before implant site preparation.
- Fluoroscopic control must be used at every step of vertebral body access.
- Exercise caution to avoid anterior wall perforation while the guidewire is inserted.

Vertebral body access

Step 1



Trocar

- Insert the trocar through the pedicle into the posterior one-third of the vertebral body.
- Remove the inner part of the trocar.
- Order trocars separately. This includes:
 - Beveled Trocar 11G
 - Diamond Trocar 11G



Guidewire

- Insert the guidewire through the trocar halfway into the vertebral body.
- Remove the tube portion of the trocar.



Tips and tricks

In wedge fractures, place implants as anterior as possible. These fractures retract slightly when opening.

Tips and tricks

There is no need to impact the guidewire into cancellous bone.

Implant site preparation

Step 2



Guidewire and reamer

- Insert the preassembled reamer/working cannula over the guidewire into the vertebra.
- Rotate the set to open the surface of the cortical bone.



Reamer

- Drill into the posterior one-third of the vertebral body.
- Remove the guidewire.
- Continue to drill until the desired position of the implant is reached.
- Disconnect the reamer from the working cannula. Unscrew and pull to remove the reamer from the vertebra, leaving the working cannula in place. The working cannula remains in place to act as a guide for the remaining instruments.




Tips and tricks

Do not change the direction of the reamer while the guidewire is still inside. This may cause the guidewire to bend and lead to difficulties removing it.

Implant site preparation

Step 2



Template

- Clean the implant site with the template.



Cannula plug

- Insert the cannula plug to the same depth as the template in order to:
 - visualise the depth of the first implant in order to position the second implant accordingly (radiopaque marker).
 - stabilise the working cannula on the first site while the second implant site is prepared.
 - stop the bleeding while the second implant site path is prepared.
- Repeat the preparation steps for the second implant site.



Tips and tricks

The instrumentation has been developed to provide full control over implant placement. The reamer and template allow for precise positioning of the implant.


Tips and tricks

Beware not to push the working cannula deeper after drilling. Next steps will be adapted to the position of the working cannula.

Implant insertion and expansion

Step 3


Prepare both implant sites before the insertion and expansion of the implants. One reamer is used to prepare both implants sites, so it should be assembled with the second working cannula. Use fluoroscopic controls at all times throughout implant site preparation.



Insert SpineJack® implants

- Insert an implant expander into each prepared path.

The orientation of the implant is gauged using the palm-held grey handle. Exercise caution to ensure that the desired orientation is achieved before expansion.




Expand SpineJack® implants

- Begin expansion on one implant by holding the palm handle with one hand and rotating the butterfly handle clockwise with the other hand.

Tips and tricks

SpineJack® system wings have been designed to allow for plastic deformations in order to adapt to the patient-specific bone conditions and vertebral endplate shape.




Reduce fracture


- Continue turning the handles until the desired vertebral body reduction is achieved.
- Use fluoroscopic controls regularly throughout the implant's expansion to ensure the desired fracture reduction is achieved.
- After each handle rotation, allow time for the bone to adjust to the implant's expansion.
- Once the expansion of the implant has started the implant cannot be closed again.

Expander removal

Step 4



- Unscrew the quick release pin counterclockwise to release the implant.



- Pull the quick release pin in the same direction, leaving the working cannula and the implant expander tube in place.

Tips and tricks

In rare cases, bone resistance may lead you to unscrew the grey handle by four or five turns before you unscrew the colored pin.

PMMA cement preparation and injection

Step 5

Use SpineJack® system in combination with the following PMMA cements:

- Cohesion® Bone Cement (CM0300)
- VertaPlex™ HV Radiopaque Bone Cement Twin Pack (0406-622-000)
- VertaPlex™ Radiopaque Bone Cement Twin Pack (0406-422-000)
- Prepare bone cement.
- TC05003/TC04003: prepare the cement with the Vexim Cement Mixing System (VCMS001).
- TC05004/TC04004: prepare the cement with:
 - Vexim Cement Injection Kit (VCIK001)
 - Masterflow™ Injection System (MF001)
 - PCD Precision System with Long 90 Extension Tube (0506-486-000)
 - Autoplex System Disposable Cement Mixer (0605-887-000)




SpineJack® cement pushers

TC04003 Cement pusher for Ø 4.2mm (6 pack)

TC05003 Cement pusher for Ø 5.0mm (6 pack)

TC05003 Cement pusher for Ø 5.8mm (6 pack)



SpineJack® injector transfer tubes

TC04004 Injector transfer tube for Ø 4.2mm (6 pack)

TC05004 Injector transfer tube for Ø 5.0mm (6 pack)

TC05004 Injector transfer tube for Ø 5.8mm (6 pack)

Internal volume

TC05003/TC05004 (5.0mm and 5.8mm) = 1.8cc

Internal volume

TC04003/TC04004 (4.2mm) = 0.9cc



Injection Systems are additional instruments and should be ordered separately.



PMMA cement preparation and injection

Step 5



Cement injection with cement pushers

- Insert the cement pusher into the working cannula/implant expander.
- TC05003/TC04003: Push the mandrel to inject PMMA cements.
- Remove empty cement pusher.
- Insert another cement pusher and continue PMMA cement injection.
- Inject PMMA cement in both sides of the vertebra.
- Continue the process until the desired quantity of PMMA cement has been injected.



Cement injection with injector transfer tube

- Insert one injector transfer tube (TC05004/TC04004) into each working cannula and clip them.
- Connect the luer lock of the chosen injection system to the injector transfer tube (TC05004/TC04004).
- Use the chosen injection system.
- When the desired quantity of PMMA cement has been injected, unclip the injector transfer tube.

Tips and tricks


Rotate the cement pushers a few times before removal. This will help prevent the creation of a “cement mouse tail” in the pedicle.

Tips and tricks

For osteoporotic VCFs, PMMA cement injection is critical to long-term results. It is recommended to inject bone cement so as to bridge the superior and inferior endplates. Use fluoroscopic controls during injection to monitor the flow.

PMMA cement preparation and injection

Step 5



- Remove the cement pusher or the injector transfer tube.
- Insert a guidewire into the hole of the upper extremity of the implant expander tube, placed inside of the working cannula.
- Simultaneously rotate and pull the guidewire and expander tube to disconnect the working cannula and expander tube from the implant.
- Remove the guidewire, working cannula and expander tube from the patient's body.



- Close the surgical access.


Clinical cases

Trauma



45-year-old patient

A.2.3 fracture in L1 after a fall from a ladder.

Hospital stay = two days



Pre-op Post-op



Pre-op and post-op superior endplate reconstructions

Trauma


57-year-old patient

A.3.2 fracture in L1 after a fall from her bed.



Low energy trauma in osteoporotic bone.

Surgery on fracture day+7.

Hospital stay = two days



Pre-op Post-op



Pre-op and post-op superior endplate reconstructions

Osteoporotic


83-year-old patient

A.3.3 fracture in L1 after a fall accident.



Traumatic VCF in osteoporotic bone.

Surgery on fracture day+10.

Hospital stay = two days



Pre-op Post-op



Pre-op and post-op superior endplate reconstructions

References

1. Oner F et al. Changes in the disc space after fractures of the thoracolumbar spine. *Journal of Bone and Joint Surgery. British volume*. 1998; 80(5):833-839.
2. Oda I et al. Does spinal kyphotic deformity influence the biomechanical characteristics of the adjacent motion segments: An in vivo animal model. *Spine*. 1999; 24(20):2139-2146.
3. Schlaich C et al. Reduced pulmonary function in patients with spinal osteoporotic fractures. *Osteoporosis International*. 1998; 8(3):261-267.
4. Lombardi I et al. Evaluation of pulmonary function and quality of life in women with osteoporosis. *Osteoporosis International*. 2005; 16(10):1247-1253.
5. Yang H et al. Changes of pulmonary function for patients with osteoporotic vertebral compression fractures after kyphoplasty. *Journal of Spinal Disorders and Techniques*. 2007; 20(3):221-225.
6. Tzermiadianos M et al. Altered disc pressure profile after an osteoporotic vertebral fracture is a risk factor for adjacent vertebral body fracture. *European Spine Journal*. 2008; 17:1522-1530.
7. Wang X et al. Kyphosis recurrence after posterior short-segment fixation in thoracolumbar burst fractures. *Journal of Neurosurgery: Spine*. 2008; 8(3):246-254.
8. Kerttula L et al. Post-traumatic findings of the spine after earlier vertebral fracture in young patients: clinical and MRI study. *Spine*. 2000; 25(9):1104-1108.
9. Cinotti G et al. Degenerative changes of porcine intervertebral disc induced by vertebral endplate injuries. *Spine*. 2005; 30(2):174-180.
10. Brinckmann P et al. The influence of vertebral body fracture, intradiscal injection, and partial discectomy on the radial bulge and height of human lumbar discs. *Spine*. 1985; 10(2):138-145.
11. Malcolm B et al. Post-traumatic kyphosis. A review of forty-eight surgically treated patients. *Journal of Bone and Joint Surgery*. 1981; 63(6):891-899.
12. Whitesides T. Traumatic kyphosis of the thoracolumbar spine. *Clinical Orthopaedics and Related Research*. 1977; 128:78-92.
13. Baeesa S et al. The efficacy of a percutaneous expandable titanium device in anatomical reduction of vertebral compression fractures of the thoracolumbar spine. *Saudi Medical Journal*. 2015; 36(1):52-60.
14. Noriega et al. Long-Term Benefits of Percutaneous Anatomical Restoration of Vertebral Compression Fractures Linked to Malignancy. *Turk Neurosurg*, 2016
15. Noriega et al. Clinical Performance and Safety of 108 SpineJack Implantations: 1-Year Results of a Prospective Multicentre Single-Arm Registry Study. Hindawi publishing corporation, 2015
16. Noriega et al. Safety and clinical performance of kyphoplasty and SpineJack® procedures in the treatment of osteoporotic vertebral compression fractures: a pilot, monocentric, investigator-initiated study. *Osteoporos Int*, 2016
17. Vanni et al. Third-generation percutaneous vertebral augmentation systems. *J Spine Surg* 2016;2(1):13-20
18. Vexim. (2005). VCF classification & guidelines for anatomical restoration [booklet].

Interventional Spine

The SpineJack® system is indicated for use in the anatomical reduction of vertebral compression fractures that may result from osteoporosis, trauma (VCF type A Magerl classification) and malignant lesions (myeloma or osteolytic metastasis).

Complications are rare. Serious adverse events, some with fatal outcome, associated with the use of bone cements for vertebroplasty, kyphoplasty and sacroplasty include myocardial infarction, cardiac arrest, cerebrovascular accident, pulmonary embolism and cardiac embolism. Although it is rare, some adverse events have been known to occur up to one year post-operatively. Additional risks exist with the use of bone cement. Please see the IFU for a complete list of potential risks.

This document is intended solely for the use of healthcare professionals.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have any questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: AutoPlex, SpineJack, Cohesion, Stryker and VertaPlex. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

2019-19721_Rev-None

Copyright © 2019 Stryker

Manufactured by:

Vexim SA,
8 rue de Vidailhan,
Hills Plaza Bât B
31130 Balma, France.

stryker.com